

NE1 Standard Analogue Water Baths

Clifton Range®



The NE1 series Clifton analogue thermostatic baths offer reliable, economical performance without compromising quality.

Comprising of a high quality corrosion resistant stainless steel tank, with a bright clean finish housed in a durable painted outer case. Painted surfaces feature anti-bacterial finish with good chemical resistance. Side lifting handles.

Features a lockable temperature dial helping prevent accidental adjustments during experiments.

	NE1-2.5	NE1-4	NE1-8	NE1-14
GENERAL				
Chamber Capacity	2.5 Litre	4 Litre	8 Litre	14 Litre
Temperature Range	Ambient + 5°C - 100°C			
Temperature Control	Analogue dial pointer, graduated in 5°C divisions			
Sensitivity/Uniformity	±0.5°C / ±0.5°C	±0.25°C / ±0.1°C		
Heater Power (220-240V)	250W	400W	800W	1000W
Voltage	230V 50-60Hz	120V or 230V models available 50-60Hz		
Drain Tap Fitted				
DIMENSIONS				
External (mm)	270w x 190d x 240h	332w x 185d x 290h	332w x 270d x 290h	361w x 332d x 290h
Internal (mm)	230w x 130d x 150h	300w x 150d x 150h	300w x 240d x 150h	325w x 300d x 150h
Water Level Max. (mm)	130			
Water Level Min. (mm)	55			

NE1-22		NE1-28		NE1-56	
GENERAL					
Chamber Capacity	22 Litre	28 Litre		56 Litre	
Temperature Range	Ambient + 5°C - 100°C				
Temperature Control	Analogue dial pointer, graduated in 5°C divisions				
Sensitivity/Uniformity	±0.25°C / ±0.1°C				
Heater Power (220-240V)	1500W			2500W	
Voltage	120V or 230V models available 50-60Hz				
Drain Tap Fitted					
DIMENSIONS					
External (mm)	538w x 332d x 290h	538w x 332d x 290h		660w x 540d x 290h	
Internal (mm)	500w x 300d x 150h	500w x 300d x 200h		620w x 500d x 200h	
Water Level Max. (mm)	130	180			
Water Level Min. (mm)	55				

NE1 Standard Analogue Water Baths

Accessories - Lids



Stainless Steel Lift Off Gable Lid

Bath	NE1-2.5	NE1-4	NE1-8	NE1-14	NE1-22	NE1-28	NE1-56
Accessory Code	N/A	SL1-4	SL1-8	SL1-14	SL1-22		SL1-56



Stainless Steel Flat Lid

Bath	NE1-2.5	NE1-4	NE1-8	NE1-14	NE1-22	NE1-28	NE1-56
Accessory Code	LD-2.5	LD-4	LD-8	LD-14	LD-22		LD-56

Please note - The NE1-56 comes complete with LD-56 2 piece lid.



Stainless Steel Hinged Gable Lid

Bath	NE1-2.5	NE1-4	NE1-8	NE1-14	NE1-22	NE1-28	NE1-56
Accessory Code	N/A			SL1-14H	SL1-22H		N/A



Clear Hinged Gable Lid

Bath	NE1-2.5	NE1-4	NE1-8	NE1-14	NE1-22	NE1-28	NE1-56
Accessory Code	N/A	GL1-4	GL1-8	GL1-14	GL1-22		N/A



Clear Lift Off Gable Lid

Bath	NE1-2.5	NE1-4	NE1-8	NE1-14	NE1-22	NE1-28	NE1-56
Accessory Code	N/A	GL1-4AF	GL1-8AF	GL1-14AF	GL1-22AF		N/A



Stainless Steel Concentric Ringed Lid

Bath	NE1-2.5	NE1-4	NE1-8	NE1-14	NE1-22	NE1-28	NE1-56
Accessory Code	N/A	SLR1-4	SLR1-8 SLR2-8	SLR1-14 SLR2-14	SLR1-22		N/A

SLR1-4: 2 x 114mm holes

SLR1-8: 4 x 83mm holes

SLR2-8: 2 x 114mm holes

SLR1-14: 4 x 114mm holes

SLR2-14: 6 x 83mm holes

SLR1-22: 6 x 114mm holes



Polypropylene Spheres (Pack of 200)

Bath	NE1-2.5	NE1-4	NE1-8	NE1-14	NE1-22	NE1-28	NE1-56
No. of Packs	1 Pack		2 Packs		4 Packs		8 Packs
Accessory Code	BP0368						

NE1 Standard Analogue Water Baths

Accessories



Stainless Steel Test Tube Racks

Bath	NE1-2.5	NE1-4	NE1-8	NE1-14	NE1-22	NE1-28	NE1-56
No. of Racks	N/A	1	2	4	6		12
Accessory Code	6870 : 26 holes x 17mm DIA						
	6871 : 16 holes x 26mm DIA						
	6872 : 36 holes x 13mm DIA						
	6873 : 18 holes x 19mm DIA						
	6900 : 12 holes x 32mm DIA						



Stainless Steel Wire Basket

Bath	NE1-2.5	NE1-4	NE1-8	NE1-14	NE1-22	NE1-28	NE1-56
Accessory Code	N/A	BSK-4	BSK-8	BSK-14	BSK-22		BSK-22*

* NE1-56 can hold 2 x BSK-22 baskets



Bath Armor™

Bath	NE1-2.5	NE1-4	NE1-8	NE1-14	NE1-22	NE1-28
No. of packs	1 x BA-5	1 x BA-5	1 x BA-12	1 x BA-12	2 x BA-12	2 x BA-12
Accessory Code	BA-5 : 4 Litre Pack					
	BA-12 : 12 Litre Pack					

Universal Accessories

CAL	Calibration Record
-----	--------------------