



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Advanced Laboratory Range

Advanced Laboratory Freezer RAFR05204

Robustly constructed with a stainless Steel exterior this freezer from the Labcold Advanced range is designed to provide close temperature control for temperature sensitive products.

Insulated with Polyurethane Water Base Insulation the refrigerant is also ammonia, HCFC and CFC free and uses an eco friendly refrigerant. Fitted with a sturdy door lock, this Advanced freezer is ideal for demanding environments such as busy laboratories and features a manual defrost putting you in control.

Remote alarm contacts and an access port are also fitted as standard.

Specifications

Product Code	RAFR05204
Gross Capacity	115 litres
Style	Undercounter
UK Warranty	2 Years
External Dimensions (mm) HxWxD	835 x 600 x 600
Net Weight	45kg
Shelves	2
Temperature Range	-15°C to -25°C

Specifications (cont.)

Factory Set Point	-20°C
Door Style	Solid
Lock Type	Key lock
Temperature Display	Yes
High/Low Temperature Alarm	Yes
Door Open Alarm	Yes
Mains Failure Alarm	No
Digital Temperature Control	Yes
Digital Temperature Display	Yes
Access Port	Yes
Access Port Size	12mm
Self Closing Door	Yes
Remote Alarm Contacts	Yes
External Material	Stainless Steel
Internal Material	Stainless Steel
Power Consumption Kwh/24hrs	1.5
Ambient Operating Temperature	16°C to 32°C

Features

- Stainless steel construction
- High temperature alarm
- Door open alarm
- Adjustable shelves
- Remote alarm contacts
- Digital temperature display
- Low temperature alarm
- Door lock
- Access port
- 2 year UK parts and labour warranty

Models

	Model	Gross Capacity	Door Style	Lock Type
	Advanced Laboratory Freezer RAFR05204	115 litres	Solid	Key lock