



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

## Sparkfree Freezers

# Sparkfree Freezer 200L RLCF0720

All Labcold chest freezers provide very stable storage conditions and slow temperature rise if power is interrupted.

They are rated for low energy consumption and include a lid lock for extra safety. Locating samples can be made simple and straight forward by use of a Labcold sample racking system – please ask us for a quote when you place your order.

## Specifications

|                                |                 |
|--------------------------------|-----------------|
| Product Code                   | RLCF0720        |
| Gross Capacity                 | 200 litres      |
| Style                          | Chest           |
| UK Warranty                    | 1 year          |
| External Dimensions (mm) HxWxD | 865 x 725 x 655 |
| Internal Dimensions (mm) HxWxD | 650 x 520 x 450 |
| Net Weight                     | 43kg            |
| Factory Set Point              | -20°C           |

## Specifications (cont.)

|                               |                  |
|-------------------------------|------------------|
| Temperature Range             | -18°C to -25°C   |
| Lid Style                     | Solid            |
| Lock Type                     | Lid lock         |
| External Material             | White Galv Steel |
| Internal Material             | Steel            |
| High Temperature Alarm        | Visual           |
| Baskets                       | 1                |
| Mains Indicator Light         | Yes              |
| Manual Defrost                | Yes              |
| Defrost Drain                 | Yes              |
| Eco Friendly Refrigerant      | Yes              |
| Power Consumption Kwh/24hrs   | 0.43             |
| Ambient Operating Temperature | 16°C to 32°C     |

## Features

- Very stable storage conditions
- Lockable
- Mains indicator light
- Operator controlled defrost
- Over temperature visual alarm (some models only)

## Models

|                                                                                    | Model                           | Gross Capacity | Lock Type |
|------------------------------------------------------------------------------------|---------------------------------|----------------|-----------|
|  | Sparkfree Freezer 200L RLCF0720 | 200 litres     | Lid lock  |
|  | Sparkfree Freezer 416L RLCF1520 | 416 litres     | Lid Lock  |
|  | Sparkfree Freezer 556L RLCF2120 | 556 litres     | Lid lock  |