



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

- 40°C Labcold Laboratory Freezers

-40°C Superfreezer 375L RLHE1345

For laboratory and general use where there is a need to store at lower than normal freezing temperatures and/or there is concern that volatile, flammable materials or components might be stored now, or at some point in the future.

Chest configuration provides improved temperature uniformity and stability with much reduced temperature rise following door openings.

Specifications

Product Code	RLHE1345
Gross Capacity	375 litres
Style	Chest
UK Warranty	1 Year
External Dimensions (mm) HxWxD	860 x 1500 x 730
Internal Dimensions (mm) HxWxD	640 x 1295 x 445
Temperature Range	-20°C to -45°C
Net Weight	82kg

Specifications (cont.)

Digital Temperature Display	Yes
Lock Type	Key lock
Lid Style	Solid
External Material	Aluminum
Internal Material	Panted Steel
High/Low Temperature Alarm	Yes
Manual Defrost	Yes
Power Consumption Kwh/24hrs	0.8
Eco Friendly Refrigerant	Yes
Ambient Operating Temperature	16°C to 32°C

Features

- Temperature stability
- CFC and HCFC free
- Operator controlled defrost
- Lockable
- Digital temperature display

Models

	Model	Gross Capacity	Lock Type
	-40°C Superfreezer 237L RLHE0845	237 litres	Key lock
	-40°C Superfreezer 314L RLHE1145	314 litres	Key lock
	-40°C Superfreezer 375L RLHE1345	375 litres	Key lock