



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Sparkfree ATEX Refrigerators and Freezers

Sparkfree ATEX Freezer RLVF1517-ATEX

This 406 litre capacity sparkfree ATEX freezer is manufactured with an interior which is licensed as a Zone hazardous area under 2014/34/EU (ATEX Directive) in order to store products which are classed in explosion group IIB+H2 and temperature class T6 under 2014/34/EU (ATEX Directive).

This free standing freezer is also fitted with a door lock, digital temperature display with alarms and backed by a free 2 year UK parts and labour warranty.

WARNING: Do not store open containers of volatile substances in this freezer

Specifications

Product Code	RLPR1517-ATEX
Gross Capacity	406 litres
Style	Free standing
UK Warranty	Two years
External Dimensions (mm) HxWxD	1875 x 600 x 700
Internal Dimensions (mm) HxWxD	1690 x 440 x 500
Temperature Range	-15°C to -25°C
Factory Set Point	-20°C

Specifications (cont.)

Internal Material	Painted Aluminum
External Material	Painted Steel
Door Style	Solid
Lock Type	Key lock
Net Weight	97Kg
Shelves	6
Temperature Display	Yes
High/Low Temperature Alarm	Yes
Door Open Alarm	Yes
Mains Failure Alarm	No
Fan Assisted Cooling	ATEX Fan
Auto Defrost	No
Manual Defrost	Yes
Rollers	4
Door Earth Connection	Yes
Condensate Collection	No
Power Consumption Kwh/24hrs	3.16
Ambient Operating Temperature	16°C to 32°C

Features

- Digital temperature display
- Fan assisted cooling for even temperature distribution
- Door lock
- 2 years UK parts and labour warranty
- High/low temperature alarm
- Door open alarm
- Door earth connection

Models

	Model	Gross Capacity	Door Style	Lock Type
	Sparkfree ATEX Fridge RLPR1517-ATEX	439 litres	Solid	Key lock
	Sparkfree ATEX Fridge RLPR1517-ATEX	439 litres	Solid	Key lock