



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

- 80°C Ultra Low Temp Freezers

-80°C ULT Freezer ULTF301

This economical and efficient chest freezer is designed to reliably store products at ultra low temperatures.

Fitted with many convenience features, such as an integral lid lock, castors and a digital temperature display with a high and low alarm, this ultra low temperature laboratory freezer even has a sparkfree interior for extra safety.

Specifications

Product Code	ULTF301
Gross Capacity	301 litres
Style	Chest
UK Warranty	1 Year
External Dimensions (mm) HxWxD	920 x 1300 x 650
Internal Dimensions (mm) HxWxD	650 x 1100 x 450
Lid Style	Solid
Net Weight	80kg

Specifications (cont.)

Digital Temperature Display	Yes
Temperature Range	-65°C to -85°C
Lock Type	Key lock
High/Low Temperature Alarm	Yes
Castors	4
Access Port	Yes
Access Port Size	12mm
Door Open Alarm	No
Remote Alarm Contacts	Yes
Sparkfree Interior	Yes
External Material	Painted Steel
Internal Material	PVC Coated Steel
Manual Defrost	Yes
Vacuum Release Port	Yes
Eco Friendly Refrigerant	Yes
Power Consumption Kwh/24hrs	7.9
Ambient Operating Temperature	16°C to 32°C

Features

- Designed for reliable temperature storage between -65 to -85oC
- Fitted with castors for easy manoeuvrability
- High and low temperature alarm
- Sparkfree interior
- Complete with access port and remote alarm contacts
- Integral lock with two keys
- Digital temperature display
- Coated steel interior
- Designed to be used with standard racks
- Environmentally friendly refrigerant

Models

	Model	Gross Capacity	Lock Type
	-80°C ULT Freezer ULTF130	130 litres	Key lock
	-80°C ULT Freezer ULTF301	301 litres	Key lock