



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Labcold Flake Ice Makers

Labcold Flake Ice Maker LITV-IQ135

With a tiny footprint compared to its ice production abilities, this ice flake maker produces soft, clean ice time after time.

Engineered from quality materials with a storage silo that is insulated to keep fresh ice in top condition this machine is ideal for a wide variety of medical applications.

Like all Labcold ice makers, it requires both a water and electricity supply and drainage 150mm below the level of the ice maker.

Specifications

Product Code	LITV-IQ135
Style	Self contained
UK Warranty	1 Year
External Dimensions (mm) HxWxD	1355 x 515 x 550
Net Weight	75kg
External Material	Stainless Steel
Ice Storage	64Kg
Daily Ice Production	164Kg

Specifications (cont.)

Automatic Production	Yes
Feet	Optional Extra
Water Use (litres per hour)	5.7
Power Consumption Kwh/24hrs	10.5
Ambient Operating Temperature	16°C to 32°C

Features

- Laboratory flake ice
 - Daily ice production 164kg
 - Stainless steel exterior
 - Air cooled
 - Eco friendly refrigerant
- Automatic production
 - Large storage bin Of up to 64kg
 - Space efficient
 - Adjustable feet are available as an optional extra

Models

Model	
	Labcold Flake Ice Maker LITV-IQ50
	Labcold Flake Ice Maker LITV- IQ85C
	Labcold Flake Ice Maker LITV-IQ135