



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Labcold Flake Ice Makers

Labcold Flake Ice Maker LITV- IQ85C

Designed for medical and laboratory purposes this low energy ice machine constructed from quality materials for the production of soft, dry flake ice.

Ideal for triage situations where ice is used to reduce swelling and pain relief, this ice maker produces clean ice flakes from drinking water. Will require access to a drain 150mm below the height of the ice maker and requires both a water and electricity supply.

Specifications

Product Code	LITV- IQ85
Style	Self contained
UK Warranty	1 Year
External Dimensions (mm) HxWxD	889 x 465 x 595
Net Weight	58kg
External Material	Stainless Steel
Ice Storage	20Kg
Daily Ice Production	88Kg

Specifications (cont.)

Automatic Production	Yes
Feet	4
Water Use (litres per hour)	3
Power Consumption Kwh/24hrs	20.1
Ambient Operating Temperature	16°C to 32°C
Eco Friendly Refrigerant	Yes

Features

- Laboratory flake ice
- 88kg daily production
- Stainless steel exterior
- Air cooled
- Eco friendly refrigerant
- Automatic production
- Large storage bin
- Space efficient
- Exterior height includes feet/legs

Models

Model

	Labcold Flake Ice Maker LITV-IQ50
	Labcold Flake Ice Maker LITV- IQ85C
	Labcold Flake Ice Maker LITV-IQ135